

FUSION FURNACE 50 PF



TECHNICAL SPECIFICATIONS

Construction

Cabinet is 5mm steel sheet mounted on a fabricated steel stand

Refractory lining

Hot face of 70% alumina firebrick backed with insulating firebrick

Internal Dimensions

400mm high x 1200mm wide x 680mm deep

Burner Equipment and Ratings

LPG – 2 x 40mm long pattern insperator 540mj/hr @ 140kpa

Natural Gas - NG2 Burner Package 540mj/hr @ 3kpa

Distillate – OIL2 Burner package 540mj/hr

Dimensions (approx)

1780 wide x 1450 deep x 1660 high

Takes 50 x 50 and 65gm, or 84 x 30,40 and 55gm Crucibles at one time.

Compatible with MARC "Multi-pour loading / unloading system

MARC provides

"Innovative and Specialised Engineered Products and Services"

Tel: 61 8 9356 4000
Fax: 61 8 9350 9333
Email: info@marcenv.com.au
Web: www.marcenv.com.au

