

Cupel Cooling Racks



Cooling Rack – 600mm

MARC manufactures cupel fuming/cooling racks in 450mm and 600mm widths for 50 and 100 cupel arrangements.

All cooling racks have 4 extracted shelves and 3 storage shelves and are made from hammertone coated mild steel.

The racks have a fume extraction spigot at the rear for direct connection to the laboratory fume extraction system.

PRODUCT DATA SHEET

Stock Number	Cupels	Dimensions (w x d x h) mm
125148	50	450 x 360 x 1408
125369	100	600 x 360 x 1408



MARC

Technologies

Ph: 61 08 9437 7400

Email: info@marctech.com.au

Website: www.marctech.com.au