

# Extracted Workstations



1200mm Workstation



Custom Double Sided Workstation

MARC manufactures a range of workstations for use in sample preparation, fire assay and metallurgical laboratories.

The benches are manufactured from mild powder coated steel tubing and heavy-duty sheet metal, with 3mm thick rubber bonded to the working surface.

Each workstation comes equipped with a fluorescent light and heavy-duty mesh shelf under the bench top illustrated.

Stock Code	Dimensions (w x d x h) mm
124986	1200 x 812 x 1652
124987	1500 x 815 x 1646
124988	1800 x 815 x 1646
127350	2000 x 815 x 1646
128775	2400 x 815 x 1646

PRODUCT DATA SHEET



**Technologies**

Ph: 61 08 9437 7400

Email: [info@marctech.com.au](mailto:info@marctech.com.au)

Website: [www.marctech.com.au](http://www.marctech.com.au)