

# Infra Red Drying Lights



Infra Red Drying Lights

MARC manufactures an eight-lamp infrared drying module suitable for the drying of metallurgical sample cakes. The module is equipped with an extraction fan to remove noxious vapors and a Perspex screen to shield the operator from the infrared lights.

Modules can be locked together to a layout of your choice, or standalone units to suit your requirements can be made.

Please note: Fan size may be effected by the length of ducting to suit the size of your shed/lab.

Stock Number	Lamps
125042	8
125043	16

PRODUCT DATA SHEET



**MARC**

**Technologies**

Ph: 61 08 9437 7400

Email: [info@marctech.com.au](mailto:info@marctech.com.au)

Website: [www.marctech.com.au](http://www.marctech.com.au)