

# Knock Up Bench

---



Typically used in fire assay laboratories, our knock up bench comes standard with an anvil for button squaring and a waste chute for easy slag disposal into a bin.

The bench has a durable Hammertone finish.



Knock Up Bench with Anvil

Stock Number	Dimensions (W x D x H) mm
126605	1258 x 760 x 895 (1045)

PRODUCT DATA SHEET



Ph: 61 08 9437 7400

Email: [info@marctech.com.au](mailto:info@marctech.com.au)

Website: [www.marctech.com.au](http://www.marctech.com.au)