

Laboratory Leach Vessels and Agitators



Agitator Stirrers

MARC manufactures Perspex laboratory leach vessels in 5L and 10L volumes for use in metallurgical laboratories.

We also manufacture a range of steel agitators to suit all needs. Special agitators can be designed and supplied on request.



Leach Vessels

Stock Number	Description
125038	5L Leach Vessel
125037	10L Leach Vessel
119230	4 Blade Rushton Radial Impeller
119229	2 Blade Radial Impeller

PRODUCT DATA SHEET



Ph: 61 08 9437 7400
Email: info@marctech.com.au
Website: www.marctech.com.au