


Mobile Benches

MARC Technologies manufactures a range of Laboratory Benches to suit many applications.



Our standard steel mobile benches are available in sizes 1200, 1500, 1800, 2000 and 2400mm. They feature a mesh shelf at half height and are finished in Hammertone enamel or 2 pack.

Depending on your requirements our mobile benches can be supplied with a variety of castors to suit your particular applications. For example preparation applications nylon/polyurethane castors provide for smooth, quiet operation and excellent wear resistance.

Benches also available in SS where required.

Non Standard sizes available on request

Stock Code	Dimensions
125247	1800mm x 900mm
129012	1500mm x 900mm

PRODUCT DATA SHEET



Ph: 61 08 9437 7400
Email: info@marctech.com.au
Website: www.marctech.com.au