

# Pneumatic Single Sieve Vibrator



Pneumatic Sieve Vibrator



We manufacture a pneumatically operated single sieve vibrator for 200mm or 450mm sieves in wet environments.

The vibrator is extremely cost effective and operates effectively.

The vibrator unit clips around a standard sieve and uses compressed air only. The vibrator may be used with water without risk of electric shock.

Please note that hearing protection as with almost all sample preparation and metallurgical equipment is strongly recommended.

The vibrator is supplied with fitting and air hose pack as shown in the inset.

Stock Number	Sieves
119 639	200mm
119 640	450mm

PRODUCT DATA SHEET



Ph: 61 08 9437 7400  
Email: [info@marctech.com.au](mailto:info@marctech.com.au)  
Website: [www.marctech.com.au](http://www.marctech.com.au)