

Pouring and Hot Pot Trolleys

PRODUCT DATA SHEET



Pouring Trolley - 50 Pot

MARC Manufactures a range of fire assay pouring and hot pot trolleys for 25, 30 and 50 pot/crucible configurations. These trolleys have been designed to clip together for simultaneous make use and are designed for use with the optional pouring tong.

Hot pot spike trays are not supplied as standard with the hot pot trolley but can be supplied on request.

The trolleys are made from mild steel and have a durable hammertone finish.



Hot Pot Trolley - 50 Pot

Pouring Trolleys	
Stock Number	Pots
125119	25
128782	30
125492	50

Hot Pot Trolleys	
Stock Number	Pots
125121	25
125090	30
125473	50

Spike Trays (to be ordered separately)	
Stock Number	Pots
125085	25
125086	30
125087	50



Ph: 61 08 9437 7400
 Email: info@marctech.com.au
 Website: www.marctech.com.au