

SS450M Sieve Shaker & Sound Enclosure



Sieve Shaker

MARC manufactures a robust sieve shaker for dry sieving with either 450mm square or round sieves.

The SS450M is a heavy-duty, floor mounted oscillating sieve shaker designed for applications where customers need a heavy-duty laboratory screen to process larger metallurgical samples.

The machine can be fitted with up to eight 450mm square or round perforated plate sieves, standard BS or ASTM perforated plates or wire mesh test sieves.

The SS450M is supplied as standard with a separate IP55 enclosed on/off control panel complete with electronic timer for installation remote from the unit.

A sound enclosure that allows access to 2 sides of the sieve shaker is a recommended accessory. If ordered together with the sieve shaker, the controls can be wired into the enclosure.



Sound Enclosure

PRODUCT DATA SHEET

Stock Number	Item	Dimensions (w x d x l) in mm
125697	SS450	900 x 1050 x 513
120617	Sound Enclosure	TBA



Ph: 61 08 9437 7400
Email: info@marctech.com.au
Website: www.marctech.com.au