


Test Tube Tumbler



Test Tube Tumbler



MARC manufactures a small bench mounted stainless steel test tube tumbler ideally suitable for small batches of samples as used in the gold industry for DIBK extractions for plant samples.

The internal holding container can be modified to take test tube racks of various sizes. Please state test tube size and quantity when requesting quote.

The unit operates on the same principle as a bottle roller and has no pinch points and passes all safety regulations.

Larger floor mounted units can also be made if required.

Stock Number	Tumbler Case Dimensions (w x h x d) in mm
128 613	410 x 204 x 214

PRODUCT DATA SHEET



MARC

Technologies

Ph: 61 08 9437 7400

Email: info@marctech.com.au

Website: www.marctech.com.au