



MARC Technologies

Nuline Refrigeration



***Health Care & Laboratory
Refrigerators/Freezers***

Laboratory Refrigerator

NLM 700/1

Standard features:

- Upright Laboratory refrigerator
- TGA Approved
- Top mounted system
- Five (5) rows of PVC coated shelves
- Interior Lighting
- Castors
- Automatic water disposal
- Digital Display with Built in Alarms
- Cabinet temperature is 2°C – 8 °C
- Operating ambient temperature is 40°C
- White colourbond interior and exterior finish with Stainless steel interior floor



Model	Dimensions w x d x h
NLM700/1	870 x 700 x 1980 plus castors / legs

Optional features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup NV Data Logger
- Voltage Free Connection
- Port entry
- Stainless Steel Interior
- Stainless Steel Exterior

Laboratory Refrigerator

NLM 1000/2, NLM1614/3 & NLM2170/4

Standard features:

- Upright Laboratory refrigerator
- Top mounted system
- TGA Approved
- Five (5) rows of PVC coated shelves behind each door
- Interior Lighting
- Castors
- Automatic water disposal
- Digital Display with Built in Alarms
- Cabinet temperature is 1°C- 12°C
- Operating ambient temperature is 40°C
- White colour bond interior and exterior cabinet finish with stainless steel interior floor.



Model	Dimensions w x d x h
NLM 1000/2	1200 x 700 x 1980 plus castors / legs
NLM 1614 /3	1880 x 700 x 1980 plus castors / legs
NLM 2170 /4	2400 x 700 x 1980 plus castors / legs

Optional features:

- Chart recorder with Battery Backup & USB Port0
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry
- Stainless Steel Interior
- Stainless Steel Exterior

NLMB Laboratory & Medical

NLMB 2

Standard Features:

- Bottom Mounted Condensing Unit
- Temperature Indicator
- Cabinet Finish in white Colourbond – Interior / Exterior with Stainless Floor
- 4 P.V.C Adjustable full width Shelves
- Double Glazed Doors
- Mullion Lighting
- Mounted on Castors
- Automatic Condensate water Disposal
- Lockable Door Handles



Model	Dimensions w x d x h	Doors
NLMB 1	750 x 700 x 1885	1
NLMB 2	1200 x 700 x 1885	2
NLMB 3	1660 x 700 x 1885	3

Optional Features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry
- Stainless Steel Interior
- Stainless Steel Exterior

NLAB – Premium

Standard Features:

- Audible Alarms
- Lockable Door Handle
- Automatic Water disposal
- Voltage Free alarm connection
- Digital Temperature indicator
- Door Ajar alar
- LED Lighting



Model	Dimensions w x d x h	Doors
NLAB 1	870 x 700 x 1980	1
NLAB 2	1200 x 700 x 1980	2
NLAB 3	1800 x 700 x 1980	3

Optional Features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry

Laboratory Refrigerator –Glass Doors



Standard features:

- Upright Display Refrigerator
- TGA Approved
- Top mounted system
- Five (5) rows of PVC coated shelves
- Castors
- Automatic water disposal
- Digital Display with Built in Alarms.
- Cabinet temperature is 2°C - 8°C
- Operating ambient temperature is 40°C
- Door Locks

Model	Dimensions w x d x h	Capacity
MF650TNG	740 x 830 (including doors) x 2010 (including castors)	650L
MF1410TNG	1480 x 830 (including doors) x 2010 (including castors)	1410L

Optional features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry

Laboratory Display Refrigerator



Standard Features:

- Upright Display Laboratory Refrigerator
- Bottom Mounted System
- 4 Rows of P.V.C Shelves
- Castors
- Automatic Defrost and Water Disposal
- Digital Display with Built in Alarms
- Cabinet Finish – White Enamelled Interior / Exterior
- Door Locks

Model	Dimensions w x d x h	Temps	Capacity
NLDR420	680x720x1990 + Castors	2°C to 8°C	412L
NLDR930	1370x720x1990 + Castors	2°C to 8°C	930L

Optional Features:

- Voltage Free Connection
- Auto Dialler with Battery Backup
- NV Data Logger
- Chart Recorder with Battery Backup & USB Port

Laboratory Freezer –Solid Doors



Standard features:

- Upright storage freezer
- TGA Approved
- Top mounted system
- Five (5) rows of PVC coated shelves
- Castors
- Automatic water disposal
- Digital Display with Built in Alarms.
- Cabinet temperature is -20°C to -25°C
- Operating ambient temperature is 40°C
- Door Locks

Model	Dimensions w x d x h	Capacity
MF650BT	740 x 830 (including doors) x 2010 (including castors)	650L
MF1410BT	1480 x 830 (including doors) x 2010 (including castors)	1410L

Optional features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry

Laboratory Display Freezer



Standard Features:

- Upright Display Laboratory Freezer
- Bottom Mounted System
- 4 Rows of P.V.C Shelves
- Interior Lighting
- Castors
- Locks
- Automatic Defrost and Water Disposal
- Digital Display with Built in Alarms
- Cabinet Finish – White Enameled Interior / Exterior
- Door

Model	Dimensions w x d x h	Temps	Capacity
NLDF420	680x720x1990 + Castors	-20°C to -24°C	412L
NLDF930	1370x720x1990 + Castors	-20°C to -24°C	930L

Optional Features:

- Voltage Free Connection
- Auto Dialer with Battery Backup
- Chart Recorder with Battery Backup & USB Port
- NV Data Logger

NHR *Spark Free* Refrigerator

Features:

- Static Refrigeration System
- Digital Display with built in Alarms
- Glass or Solid Door
- Auto Water Disposal
- Castors
- White Colourbond Interior / Exterior
- Door Lock



Model	Dimensions w x d x h	Number of Shelves	Capacity
NHRS 400	600 x 620 x 1855 + Castors	5	350L
NHRS 600	775 x 750 x 1900 + Castors	5	590L

Optional Features:

- Voltage Free Connection
- Auto Dialler with Battery Backup
- Stainless Steel Interior
- Stainless Steel Exterior
- Stainless Steel Shelves

NLM Laboratory Refrigerators – *Spark Free*

Features:

- Static Refrigeration System
- Digital Display with Built in alarms
- Auto Defrost
- Castors
- Glass or Solid Doors



Model	Dimensions w x d x h
NLMS 700/1	870 x 700 x 1980 plus castors / legs
NLMS 1000/2	1200 x 700 x 1980 plus castors / legs
NLMS 1614/3	1880 x 700 x 1980 plus castors / legs
NLMS 2170/4	2400 x 700 x 1980 plus castors / legs

Optional Features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry
- Stainless Steel Interior
- Stainless Steel Exterior

NLMBS – *Spark Free* Refrigerator

Features:

- Static Refrigeration System
- Digital Display with Built in alarms
- Auto Defrost
- Castors
- Glass or Solid Doors



Model	Dimensions w x d x h	Doors
NLMBS 2	1200 x 650 (plus doors) x 1900	2
NLMBS 3	1600 x 650 (plus doors) x 1900	3

Optional Features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialer with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry
- Stainless Steel Interior
- Stainless Steel Exterior

Spark Safe Freezer

HF Series



HF200



HF400



HF600

Features:

- Digital Display with Built in Alarms – High and low read out
- Static Refrigerated Shelves
- Door Lock
- Bottom Mounted System
- TGA Approved

Model	Capacity	Dimensions w x d x h	Number of Shelves
HF200	135 Lt	600x620x850	3
HF400	350 Lt	600x620x1900	5
HF600	570Lt	770x750x1900	5

Optional Features:

- Voltage Free Connection
- Data Logger NV600
- Auto Dialler
- Port Entry

NHF – Spark Safe Freezer

Features:

- Static Refrigeration System
- Digital Display with Built in alarms
- Auto Defrost
- Castors
- Stainless Steel Interior



Model	Capacity	Dimensions w x d x h	Number of Shelves
NHF400	356L	650 x 750 x 1855*	5
NHF600	505L	800 x 750 x 1900*	5

*Height is plus castors or legs

Optional Features:

- Voltage Free Connection
- Data Logger NV600
- Auto Dialler
- Port Entry

NHRTS – Top Mounted *Spark Safe* Refrigerator

Features:

- Static Refrigeration System
- Digital display with built in alarms
- Auto defrost
- Castors
- Top mounted system
- Stainless steel interior



Model	Capacity	Dimensions w x d x h	Number of Shelves
NHRTS 400	350L	600 x 620 x 1855*	5
NHRTS 600	570L	775 x 750 x 1900*	5

*Height is plus castors or legs

Optional Features:

- Voltage Free Connection
- Data Logger NV600
- Auto Dialler
- Port Entry

NLMTS – Top Mounted *Spark Safe* Refrigerator

Features:

- Static Refrigeration System
- Digital display with built in alarms
- Glass or Solid Doors
- Castors
- Auto Defrost



Model	Dimensions w x d x h
NLMTS 700/1	870 x 700 x 1980 plus castors / legs
NLMTS 1000/2	1200 x 700 x 1980 plus castors / legs
NLMTS 1614/3	1880 x 700 x 1980 plus castors / legs
NLMTS 2170/4	2400 x 700 x 1980 plus castors / legs

Optional Features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry
- Stainless Steel Interior
- Stainless Steel Exterior

2 Temp Refrigerator / Freezer



Features:

- Digital Display with Built in Alarms
- Auto Defrost in 2°C Section
- Manual Defrost in -20°C Section
- Stainless Steel Interior / Exterior
- Top Mounted System

Model	Dimensions w x d x h	Capacity 2°C/ Shelves	Capacity -20°C / Shelves
NDT	870 x 750 x 2120 (including castors)	386L / 3 Shelves	212L / 2 Shelves

Optional Features:

- Auto Dialler with Battery Backup
- Voltage Free Connection
- NV Data Logger
- Stainless Steel Shelves

HRF Combination

HRF 400 2T

Features

- QCPP Approved
- Digital Temperature indicator
- Built in Hi/Low alarms
- Door locks
- Data Logger Included – For refrigerator section only



Model	Capacity	Dimensions w x d x h	Shelves
HRF 400 2T	Refrigerator: 135L Freezer: 135L	600 x 620 x 1755 (on castors)	Refrigerator: 4 Adjustable Freezer: 3 Fixed Static

Description Cost

- Voltage Free Connection
- Data Logger NV600
- Auto Dialler
- Port Entry

HR Pharmacy Refrigerator Solid Door



HR200



HR400



HR600

Features:

- Built In Alarm
- QCPP Approved
- TGA Approved
- Digital Temperature Indicator
- Fan Forced
- Door Lock
- Data Logger Included

Model	Dimensions	Door	Shelves
HR200	600W x 600D x 850H	Solid	4
HR400	600W x 600D x 1850H	Solid	6
HR600	780W x 690D x 1900H	Solid	6

Optional Features:

- Voltage Free Connection
- Auto Dialler
- Port Entry

HR Pharmacy Refrigerator Glass Door



HR200G



HR400G



HR600G

Features:

- Built In Alarm
- QCPP Approved
- TGA Approved
- Digital Temperature Indicator
- Fan Forced
- Door Lock
- Interior Lighting
- Data Logger Included

Model	Dimensions	Door	Shelves
HR200G	600W x 600D x 850H	Glass	4
HR400G	600W x 600D x 1850H	Glass	6
HR600G	780W x 690D x 1900H	Glass	6

Optional Features:

- Voltage Free Connection
- Auto Dialler
- Port Entry

NHR – Pharmacy Refrigerator

Features:

- Built in Alarm
- Digital Temperature Indicator
- Fan Forced
- Door Lock
- Interior Lighting
- Door Lock on Handle



Model	Dimensions	Capacity	Shelves
NHR400	600 x 620 x 1855 (not including castors)	350L	6
NHR600	775 x 750 x 1900 (not including castors)	590L	6

Optional Features:

- Voltage Free Connection
- Data Logger NV600
- Auto Dialler
- Port Entry

NLM Pharmacy

Standard features:

- Upright Laboratory refrigerator
- Top mounted system
- Five (5) rows of PVC coated shelves behind each door
- Castors
- Automatic water disposal
- Digital Display with Built in Alarms
- Cabinet temperature is 2°C- 8°C
- Operating ambient temperature is 40°C
- White colourbond interior and exterior cabinet finish with stainless steel interior floor.
- Door Locks



Model	Dimensions
NLMP 700/1	870 x 700 x 1980 plus castors / legs
NLMP 1000/2	1200 x 700 x 1980 plus castors / legs
NLMP 1614 /3	1880 x 700 x 1980 plus castors / legs
NLMP 2170 /4	2400 x 700 x 1980 plus castors / legs

Optional features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry
- Stainless Steel Interior
- Stainless Steel Exterior

NLMB Pharmacy

Standard Features:

- Bottom Mounted Condensing Unit
- Temperature Indicator
- Cabinet Finish in white Colourbond – Interior / Exterior with Stainless Floor
- 4 P.V.C Adjustable full width Shelves
- Double Glazed Doors
- Mullion Lighting
- Mounted on Castors
- Automatic Condensate water Disposal
- Full Length Door Handles
- Door Locks



Model	Dimensions w x d x h	Doors
NLMBP 1	750 x 700 x 1885	1
NLMBP 2	1200 x 700 x 1885	2
NLMBP 3	1660 x 700 x 1885	3

Optional Features:

- Chart recorder with Battery Backup & USB Port
- Phasefale Temp with Battery Backup – Dual Temp
- Auto Dialler with Battery Backup
- NV Data Logger
- Voltage Free Connection
- Port entry
- Stainless Steel Interior
- Stainless Steel Exterior